



# **Product Information : StemFit™ Basic03**

#### 1. Introduction

StemFit™ Basic03 is an animal-origin free, defined medium for human pluripotent stem cell culture. This medium is capable of maintaining human ES and iPS cells under feeder-free culture conditions.

## 2. Materials Provided 1

Component	Volume	Storage
Liquid A	400ml	Store at below -20 °C
Liquid B	100ml	Store at below -20 °C



### 3. Media Preparation

StemFit™ Basic03 is provided frozen as a 2 component set containing "Liquid A" and "Liquid B" and can be stored at below -20 °C until use. Use sterile techniques to prepare StemFit™ Basic03 medium.

- 1) Before use, thaw the frozen "Liquid A" and "Liquid B" with occasionally mixing at room temperature (15-25 °C). Do not thaw, "Liquid B" at 37 °C, as it accelerates the degradation of medium ingredients. 2
- 2) Aseptically mix medium components by adding the full volume of "Liquid B" to "Liquid A". Mix thoroughly. If prepared aseptically, StemFit™ Basic03 is ready for use. If precipitations are observed, keep the bottle at room temperature and dissolve them.
- 3) Upon thawing, StemFit™ Basic03 medium may be aseptically aliquoted and stored at below -20 °C. 3 4
- 4) Before use, warm aliquots to room temperature and use immediately.



This product does not contain basic fibroblast growth factor (bFGF)



Do not thaw, "Liquid B" at 37 °C, as it accelerates the degradation of medium ingredients.



Thawed StemFit<sup>™</sup> Basic03 medium may be stored at 2-8 °C for up to two weeks.



We recommend storing the medium protected from light

#### 4. Precaution and Disclaimer

StemFit™ BasicO3 is for research use only and is not for diagnostic use, therapeutic use, or as a food.

# 5. Contact the following department for product information:

Amino Acids Dept., AminoScience Division, AJINOMOTO CO., INC.

1-15-1 Kyobashi, Chuo-ku, Tokyo 104-8315, Japan

TEL: +81-3-5250-5040 E-mail: **stemfit@asv.ajinomoto.com**