



## ● Cell expansion after thawing



Cell stocks of 1210B2 hiPSC were prepared with StemFit (red bar), Medium A (blue bar), or Medium B (green bar). Cell stocks were thawed with each respective medium and cultured for a week. Cell stocks prepared with StemFit showed superior cell growth compared to Medium A or Medium B. The superior cell growth of StemFit allows the expansion of the cell stocks with greater efficiency.

## ● Karyotype analysis after 52 passages



201B7 cells were passaged in StemFit medium for 52 passages. G-band staining showed cells maintained a normal karyotype after 52 passages.

Eat Well, Live Well.



For further information, please contact


✉ [stemfit@asv.ajinomoto.com](mailto:stemfit@asv.ajinomoto.com)

AJINOMOTO CO., INC. AminoScience Division

15-1, Kyobashi 1-Chome, Chuo-Ku, Tokyo 104-8315, Japan



## Feeder-free medium for ES/iPS cells



# StemFit Technical tips-3

Key Points for *making cell stocks*  
with *single-cell passaging method*



**Highly efficient cell stock production with a single-cell passaging protocol**




**Superior cell growth after thawing**


Eat Well, Live Well.



## Highly efficient cell stock production



## Excellent cell growth and stable karyotype



## Efficient protocol for making cell stocks

### Freezing

#### Freezing Protocol

##### Procedure (6-well plate)

###### < Cell preparation >

- 1 Culture cells with StemFit for more than 2 passages.
- 2 Use 50-80% confluent cells for cell stocks.

###### < Procedure >

- 1 Detach and count cells.
- 2 Aspirate the medium and resuspended the cells at  $1 \times 10^6$  cells/mL in STEM-CELLBANKER® freeze solution. (Bambanker® can be used as well)
- 3 Dispense 200  $\mu$ L ( $2 \times 10^5$  cells) per cryovial.
- 4 Cool using a programmed freezer ( $-1^\circ\text{C}/\text{min}$ ), or in a Freezing Container (e.g. CoolCell® or Mr. Frosty™) at  $-80^\circ\text{C}$  for at least 3 hrs.
- 5 Transfer to liquid nitrogen within a few days for extended storage.


Aspirate the medium and wash once with 2 mL of PBS




Add 500  $\mu$ l/well of Accutase and incubate at  $37^\circ\text{C}$  for 10 min

\* TrypLE™ can also be used for cell dissociation  
\* Incubation time may vary depending on the matrix

• Before incubation with Accutase



• Gaps in the colonies appear and dissociation of colonies is apparent



Whole colony can be smoothly detached

Count the cells and centrifuge the tubes (300 X g) at room temperature



### Thawing

#### Thawing Protocol

##### Procedure (6-well plate)

- 1 Thaw frozen stock immediately in a water bath at  $37^\circ\text{C}$  for  $\sim 1$  min until only small ice particles remain.
- 2 Transfer the cell suspension to 5 mL of culture media with gentle pipetting (1-2 times).
- 3 Spin the cells at 300 x g, RT, for 4 min.
- 4 Resuspend in 1 mL of medium and count cells.
- 5 Seed 65,000 live cells to 1.5ml culture media containing  $10 \mu\text{M}$  Y-27632.
- 6 After  $> 24$ hr, change to regular culture medium (NOT containing Y-27632).